


MICRO MIX

Sample Mixer



Main Features :

- The mixer used for mixing sample tubes, EDTA tubes, ESR Tubes.
- The continuous speed control ensures an optimum blending of the sample material.
- The sample plate can be equipped with a holding ring for further sample tubes with a different diameter.

Specifications :

| | | |
|----------------|---|--------------------------|
| Connected load | : | 110 / 220 V x 50 - 60 Hz |
| Speed | : | 10 rpm |
| Capacity | : | 16 tubes. |
| Dimensions | : | 25 x 25 x 15 cms. |
| Weight | : | 1 kg. |

MICRO LAB 
INSTRUMENTS