

WRX - 806

Elisa Plate Washer



Main Features :

- Flexibility and precision in microplate strip washing.
- Top loading : Provides easy handling and operation.
- Low noise pump.
- Optional liquid level sensors.
- Aspiration time : Continuous aspiration reduces risks of well to well overflow and cross contamination.
- Compact Size.

Specifications :

- Manual Microplate washer with 8 position manifold.
- SS Probes for dispensing buffer & aspirating wash solution.
- Manifold : 1x8 Channel.
- Bottles : Buffer 1000 ml., Waste : 1000ml.
- Residual Volume : 2 μ l (typical)
- CV : <3%

Power : 110 - 220V, 50 - 60 Hz

Operating Temp : 10° to 40°c

Dimensions : 20 x 26 x 26 cms

Weight : 4 kgs.